

What's **new** from NSK?

NEW PRODUCTS

NSK has what you need **n_xt**
next smile next achievement next innovation



next treatment



next surgery



next sterilization



next solution



MORE INFO

ORAL HYGIENE



Varios Combi Pro2

ULTRASONIC SCALER X POWDER JET COMBI DEVICE

We've designed the **Varios Combi Pro2** to deliver the advanced functionality dental professionals truly need.

Combining subgingival and supragingival powder therapy with ultrasonic scaling, it adapts to your preferred treatment style for **seamless periodontal care**.

Backed by NSK's expertise, the Varios Combi Pro2 is built for durability, effortless maintenance, and uninterrupted performance—your perfect 'daily driver' for smooth, efficient treatment.

With NSK, you have the tools to perform at your best—delivering exceptional care, from prevention to treatment, every time.

Key Features:

- **Jet nozzle** for precise powder therapy
- **Water heating** for enhanced comfort
- **Wireless foot-control** for hands-free operation
- **Titanium ultrasonic scaler handpiece** for durability and precision

INSTRUMENTS



Ti-Max Z99L

1:5 CONTRA-ANGLE

NSK's flagship model with uncompromising quality, designed for the high demands of dental professionals. With an improved chuck and gear that make this handpiece twice more durable, enhanced cutting performance, as well as meticulous refinements to the head angle, water spray angle, and bur retention force, the Ti-Max Z99L delivers an unrivaled treatment experience.



S-Max M94L

1:4 CONTRA-ANGLE

Advanced model specially designed for harder material cutting typified by Zirconia. Zirconia is an ideal material for dental restorations due to its strength and aesthetics, but its high hardness makes removal time-consuming and challenging. New 1:4 speed increasing contra-angle delivers high torque for smooth and efficient removal of zirconia.

STERILIZATION



iClave 24

CLASS B AUTOCLAVE

The new iClave 24 autoclave perfectly combines performance, user experience, and technology, ensuring an effective, comfortable, and worry-free experience. Its 24-litre stainless steel chamber with 5 trays of 2 different dimensions maximises the actual capacity. In addition, the new high-efficiency steam generator allows for superior performance, resulting in perfectly sterilised equipment in less time.


SURGICAL



SG12

ULTRASONIC BONE SURGERY TIP

The long tip design gives easy access to the mandibular ramus. Blades on all sides allow access from any angle, with three-dimensional cutting of the blades edge, allowing high cutting efficiency. The new SG12 tip inherits the cutting power of the well renowned SG1, along with greater accessibility and convenience for general bone surgery.



M-OC1L, M-OC1R, M-RE2, M-ST1

ULTRASONIC BONE SURGERY 'MAXI REACH' TIPS

These tips, approximately twice as long as existing ultrasonic bone surgery tips, allowing access to deep areas, previously inaccessible. Also its low tip height makes cutting possible even in constricted areas.